2016 2017



Dental Lamp Catalogue

UE United Consolidated Sdn. Bhd.

(Co. No. 815185H)

GST No. 001400324096

11-1, Ground Floor, Jalan PJU 1/42A,

Dataran Prima (Block E2),

47301 Petaling Jaya,

Selangor DE, Malaysia

Tel: 03 7804 6908

Fax: 03 7804 6480

Email: sales@ueunited.com.my

Web: www.ueunited.com.my

Dental Lamps

| Lamp | Code | Volts & Watts | Description |
|------|---------------------|---------------------|---|
| | P3/10 (13165) | 14V 35W | Philips MR11 halogen lamp with reflector and GZ4 base for Sirona C7 / C8 |
| 4 | 13298 | 10V 52W | Philips MR11 halogen lamp with reflector and GZ4 base |
| | 14623 | 17V 95W | Philips Focusline halogen capsule with G6.35 base. |
| | 13865 | 12V 75W | Halogen bulb with reflector G5.3 base for Lightcure and Spectrum Curing Light |
| | HLX64617 (M117) | 12V 75W | Halogen bulb with reflector G5.3 base for dental curing |
| | A1/51 | 12V 100W | Halogen bulb with reflector G5.3 base for dental curing |
| | HLX 64642 (M184) | 24V 150W | Halogen capsule bulb G6.35 base |
| | A453 | 12V 55W | Halogen H3 PK22S |
| | A1/216 | 24V 150W | Halogen capsule with G6.35 base |
| | A1/220 | 12V 50W | Halogen capsule with G6.35 base |
| | A1/230 | 12V 75W | Halogen lamp with reflector GZ6.35 base. |

Due to vast variety of lamps and light source use in medical field, we are unable to list out all that is related.

However, if there are any lamps that are not found in this catalogue, please feel free to contact us.

CALL : 03 7804 6908

FAX : 03 7804 6480

WHATSAPP : 012 3735153 - MR. ALEX

WHATSAPP : 012 2556479 - MS. TAN

EMAIL : sales@ ueunited.com.my

Provide us any of the followings, we **WILL** find the lamps for you.

- PICTURE WITH DIMENSION
- LAMP WATTAGE & VOLTAGE
- EQUIPMENT'S MADE & MODEL

While every efforts is being made to ensure that the information is as accurate as possible, as the products, photos, dimensions, diagram and data in this catalogue have been provided in good faith with help and support by our suppliers. However, we shall not be held responsible for any unforeseen errors that happen after the print.

| <u>Note</u> | |
|-------------|--|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |